

**Title: FLUOROMETRIC DETERMINATION OF CHLOROPHYLL A**

---

Author: \_\_\_\_\_ Date: \_\_\_\_\_  
Marie E. DeLorenzo

Program Manager: \_\_\_\_\_ Date: \_\_\_\_\_  
Michael H. Fulton

Branch Chief: \_\_\_\_\_ Date: \_\_\_\_\_  
Geoffrey I. Scott

---

**1.0 OBJECTIVE**

Chlorophyll *a* is used to estimate phototrophic biomass. The purpose of this method is to quantify chlorophyll *a* concentration from water samples. This method was adapted from Glover and Morris (1979).

**2.0 HEALTH AND SAFETY**

Personnel should wear lab coats and chemical resistant gloves.

**3.0 PERSONNEL/TRAINING/RESPONSIBILITIES**

Personnel should not perform this method until training by experienced individuals is complete.

**4.0 REQUIRED AND RECOMMENDED MATERIALS**

20 mL plastic scintillation vials  
acetone  
magnesium carbonate ( $\text{MgCO}_3$ )  
deionized water  
glass fiber filters (Type GF/F, 25 mm diameter)  
filter apparatus  
filter forceps  
fluorometer (e.g. Sequoia-Turner Model 450)

disposable, borosilicate glass test tubes for fluorometer

Vortex mixer

1 and 10 mL pipettes/bulbs

## **5.0 PROCEDURE**

### **5.1 Chlorophyll extraction**

- \$ Sample volume required will vary, depending on the chlorophyll concentration in the water. For PFU samples, 10 mL is typically adequate.
- \$ Filter water sample onto glass fiber filter (Type GF/F, 25 mm diameter).
- \$ Just before all sample passes the filter, rinse the column with two separate aliquots (1.0 mL each) deionized water. Continue vacuum until all liquid is gone.
- \$ Release vacuum, disassemble filter tower apparatus, and remove filter with forceps.
- \$ Place filter face up on bottom of scintillation vial, add 1 mL  $\text{MgCO}_3$  and freeze until analysis.
- \$ Samples should be kept in the dark for the rest of the procedure. To extract the samples, add 9 mL of acetone to each scintillation vial and shake well.
- \$ Refrigerate samples overnight in the dark at 4 EC, shake the samples the next day, and refrigerate overnight again.
- \$ The next day, bring the samples to room temperature and read on a fluorometer.

### **5.2 Fluorometric measurement**

#### 5.2.1 Chlorophyll *a*

- \$ Before using the fluorometer for unknowns, a standard curve should be created with pure chlorophyll *a* extracts (available from Sigma).
- \$ AZero@each door opening of the fluorometer immediately prior to use, using a tube of 90% acetone.

Effective Date: 06 December, 2000

- \$ Decant the chlorophyll extract from the scintillation vial into the fluorometer test tube. Care should be taken not to transfer particulates from the filter into the tube.
- \$ Wipe off the sides of the test tube with a Kim Wipe and place the test tube in the fluorometer.
- \$ Record the fluorescence units, door setting and gain setting.
- \$ Samples that are too concentrated may be diluted with 90% acetone.
- \$ A new test tube should be used for each sample.

#### 5.2.2 Phaeo-pigments

This extra step is performed to determine and correct for the concentration of phaeo-pigments (chlorophyll degradation products) in the samples. It is usually not necessary with the GF/F filtering technique, but is provided here for those using other collection methods.

- \$ After the first reading is taken on the fluorometer, remove the tube and add 2 drops of 5% v/v hydrochloric acid.
- \$ Mix contents of tube with a vortex mixer.
- \$ Take a second reading 30-60 seconds later, after a stable value is reached.

#### **5.3 Calculations**

Chlorophyll *a* (µg/mL)= (door factor\*(chlorophyll fluorescence reading-phaeo-pigment reading)\*gain correction\*acetone volume)/volume filtered\*gain.

### **6.0 QUALITY CONTROL/QUALITY ASSURANCE**

A minimum of three replicates per site or treatment is recommended. It is important that each sample is well-mixed prior to filtration and that the samples are kept in the dark after collection onto the filters.

### **7.0 REFERENCES**

Glover H.E. and Morris I. 1979. Photosynthetic carboxylating enzymes in marine phytoplankton. *Limnol Oceanogr* 23:510-519