

Title: **Protein Assay Protocol**

Author:	_____	Date: _____
	Pete Key	
Program Manager:	_____	Date: _____
	Michael H. Fulton	
Branch Chief:	_____	Date: _____
	Geoffrey I. Scott	

1.0 OBJECTIVE

To conduct protein assay.

2.0 HEALTH AND SAFETY

Personnel should wear a lab coat and chemical resistant gloves. Personnel should be aware that a phenolic compound is used in this assay and due caution should be taken.

3.0 PERSONNEL/TRAINING/RESPONSIBILITIES

Any employee who routinely works in the laboratory should be capable of performing this task. Training of new staff should be carried out under supervision of an experienced technical employee familiar with this SOP before the employee can work unsupervised.

4.0 REQUIRED AND RECOMMENDED MATERIALS

This section lists the required supplies and equipment:

spectrophotometer	cuvettes
assay reagents, buffers, standard	pipetors
lab coat	test tubes
chemical resistant gloves	

5.0 PROCEDURE

5.1 Protein Assay

1. Make Lowry Reagent, Folin & Ciocalteu's Phenol Reagent and Protein Standard according to **Protein Assay Notebook** in Room 230.
2. Prepare standards and blank for 50 μ l and 100 μ l protein samples according to Table in Protein Assay Notebook in Room 230.
3. Prepare Samples {See **Protein Assay Notebook** in Room 230}
 - a.) 50 μ l samples: Dilute samples with 950 μ l D.I. water to acquire 1.0 ml total sample volume
100 μ l samples: Dilute samples with 900 μ l D.I. water to acquire 1.0 ml total sample volume
 - b.) Add 1.0 ml Lowry Reagent to ALL tubes (Blanks, Standards, Samples).
Do not start mixing until **last** tube is filled.
Then vortex, then let stand for 20 minutes at room temperature.
 - c.) Add 500 μ l of Folin & Ciocalteu's Phenol Reagent to ALL test tubes and vortex after each addition.
Start 30 minute countdown for each one after each has been mixed.
 - d.) Transfer solutions to cuvettes and read in spectrophotometer at 500 nm.
4. Discard all reagents, standards and samples with plenty of water after use.

5.2 Spectrophotometer Setup for Protein Assay

1. Turn on Spec unit, then module, then printer
2. Key in "500", Press "2nd" then "Go to λ "
3. Insert Blank - Pour blank into cuvette and place into spectrophotometer
4. Press "2nd" then "Zero"
5. After zeroing, set for Standard Mode
Press "5" then "2nd" then "Select"
SAVE TESTS "no"
STDS "yes"
THRU 0 "no"
STDS Enter number of standards then Press "yes"
6. Insert Standards. Remove blank and pour standard into cuvette
7. Press "yes" or "send"
8. Enter concentration of Standard then "yes" or "send".
Repeat steps 6 and 7 until all Standards have been entered.
- a. Check printout for slope value. Should 1.00 or very close to 1.00.

Keep blank and standard cuvettes in order.

9. ACCEPT NAME Press "yes" to accept
 Press "no" to change letters.
 Use "yes" and "no" to move between letters.
 Press "Send" to choose letters; "Stop" when finished.
10. ID # Press "yes" (to change ID see step 9)
11. WV=500 "yes"
 INIT DELAY "0" then "Send"
 LO "yes"
 HI "yes"
 SLOPE "yes"
 INTERCEPT "yes"
 SIG AVG "yes"
12. Insert each sample into spec chamber and press "Run" on the module.

13. Check protein readout on printout. If numbers are greater than 10% or less than 10% of the two standards used, repeat steps 3 – 12 with standards that bracket the sample readout.

6.0 QUALITY CONTROL/QUALITY ASSURANCE

Personnel should adhere to good laboratory practices performing this assay. This procedure should always be performed with proper precautions to minimize personnel exposure to reagents.

7.0 REFERENCES

Lethal and sublethal effects of malathion on three life stages of the grass shrimp, *Palaemonetes pugio*. P. Key, M. Fulton, G. Scott, S. Layman and E. Wirth. 1998. Aquatic Toxicology. 40:311-322.

Sigma Chemical Company. 1989. Sigma diagnostics protein assay (Procedure No. P 5656). St. Louis, MO: Sigma Chemical Co.

DeWoskin, R.S. 1984. Good laboratory practice regulations: a comparison. Research Triangle Institute, Research Triangle Park, North Carolina. 63 pp.

USEPA. 1979. Good laboratory practice standards for health effects. Part 772 - Standards for development of test data. Fed. Reg. 44:27362-27375, May 9, 1979.

USEPA. 1980. Physical, chemical, persistence, and ecological effects testing; good laboratory practice standards (proposed rule). 40 CFR 772, Fed. Reg. 45:77353-77365. November 21, 1980.